
Product Datasheet

**Human MICB Protein, His Tag (MALS verified)
(orb1737427)**

Description

Human MICB Protein, His Tag is

(HEK293). It contains AA Ala 23 - Thr
...**Reactivity**

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

MICB

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 34.5 kDa. The protein migrates as 45-55 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human MICB Protein, His Tag (orb1737427) is expressed from human 293 cells (HEK293). It contains AA Ala 23 - Thr 308 (Accession # NP_005922.2).

MW

34.5 kDa

Uniprot ID

NP_005922.2

Expiration Date

6 months from date of receipt.