
Product Datasheet

IRAK1BP1 antibody [N1C3] (orb1741778)

Description

IRAK1BP1 antibody [N1C3]

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

Recombinant protein encompassing a sequence within the center region of human IRAK1BP1. The exact sequence is proprietary.

Target

interleukin 1 receptor associated kinase 1 binding protein 1

Preservatives

0.01% Thimerosal. 0.1M Tris-Glycine (pH7), 10% Glycerol

Form/Appearance

Liquid

Concentration

1 mg/ml (Please refer to the vial label for the specific concentration.)

Storage

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

29

Uniprot ID
[Q5VVH5](#)
Dilution Range

IHC-P: 1:100-1:1000. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.

Expiration Date

12 months from date of receipt.