

Product Datasheet

Taq Plus DNA polymerase (orb1742794)

Description	Taq Plus DNA polymerase
Species/Host	E. coli
Preservatives	10 mM Tris-HCl, pH 7.4, 100 mM KCl, 1 mM DTT, 0.1 mM EDTA, 0.5% Tween 20, 0.5% IGEPAL CA-630, 50% glycerol
Form/Appearance	liquid
Storage	Store product at -30°C to -15°C.
Note	For research use only
Source	Thermus aquaticus
Expiration Date	6 months from date of receipt.