

Product Datasheet

Super Mix For qPCR (+gDNA Wiper) (orb1742795)

Description	Super Mix For qPCR (+gDNA Wiper)
Form/Appearance	liquid
Storage	Store product at -30°C to -15°C.
Note	For research use only
Source	Taq Plus DNA polymerase consists of a mixture of Taq DNA polymerase and an enzyme containing 3'→5' exonuclease activity. The fidelity of Taq Plus is 6X greater than that of Taq DNA Polymerase.
Expiration Date	6 months from date of receipt.