

## Product Datasheet

# ATP-dependent DNA helicase RecQ, Active (orb1742810)

<b>Description</b>	ATP-dependent DNA helicase RecQ, Active
<b>Preservatives</b>	Recombinant protein stored in 50mM sodium phosphate, pH 8.0, 300mM NaCl, 150mM imidazole, 1 mM DTT, 10% glycerol.. 50 mM sodium phosphate, pH 8.0, 300mM NaCl, 150 mM imidazole, 1mM DTT, 10% glycerol
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	0.5 µg/µl
<b>Storage</b>	Store at -70°C
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Source</b>	E. coli
<b>MW</b>	Approx. 70kDa
<b>Expiration Date</b>	6 months from date of receipt.