

## Product Datasheet

### Proteinase K, Active (MBG) (orb1742834)

<b>Description</b>	Proteinase K, Active (MBG)
<b>Preservatives</b>	50mM Tris, 2mM MgCl <sub>2</sub> , 20mM NaCl (pH 8.0), and 50% Glycerol.. 50 mM Tris-HCL, pH 7.8, 3mM CaCl <sub>2</sub> , 50% glycerol
<b>Form/Appearance</b>	lophalized powder
<b>Storage</b>	Store product at -20°C
<b>Note</b>	For research use only
<b>Source</b>	Pichia pastoris
<b>MW</b>	29 kDa
<b>Expiration Date</b>	6 months from date of receipt.