

Product Datasheet

TtCsm6, Active (orb1742844)

Description	TtCsm6, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Thermus thermophilus
MW	52.7
Uniprot ID	Q53W17
Expiration Date	6 months from date of receipt.