

Product Datasheet

TccCas13a, Active (orb1742845)

Description	TccCas13a, Active
Preservatives	Recombinant protein stored in 150mM imidazole, 50mM sodium phosphate pH 7.5, 10% glycerol, 1mM DTT.. Recombinant protein stored in 150mM imidazole, 50mM sodium phosphate pH 7.5, 10% glycerol, 1mM DTT.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Thermoclostridium caenicola
MW	144.7
Expiration Date	6 months from date of receipt.