

## Product Datasheet

### SpCas9-NLS, Active (orb1742847)

<b>Description</b>	SpCas9-NLS, Active
<b>Preservatives</b>	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	0.1 µg/µl
<b>Storage</b>	Store at -70°C.
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Source</b>	Streptococcus pyogenes
<b>MW</b>	161.0
<b>Uniprot ID</b>	<a href="#">Q99ZW2</a>
<b>Expiration Date</b>	6 months from date of receipt.