

Product Datasheet

SpCas9, Active (orb1742848)

Description	SpCas9, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 8.0, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 8.0, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Streptococcus pyogenes
MW	161.0
Uniprot ID	Q99ZW2
Expiration Date	6 months from date of receipt.