

Product Datasheet

LwCas13a, Active (orb1742853)

Description	LwCas13a, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 8.0, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 8.0, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Leptotrichia wadei
MW	139.8
Uniprot ID	U2PWF1
Expiration Date	6 months from date of receipt.