

Product Datasheet

LbuCas13a, Active (orb1742854)

Description	LbuCas13a, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 1mM DTT, 10% glycerol.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	tag-free
Note	For research use only
Source	Leptotrichia buccalis
MW	138.5
Uniprot ID	C7NBY4
Expiration Date	6 months from date of receipt.