

Product Datasheet

EiCsm6, Active (orb1742855)

Description	EiCsm6, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 200mM imidazole, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 200mM imidazole, 1mM DTT, 10% glycerol.
Form/Appearance	liquid
Concentration	0.5 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Enterococcus italicus
MW	51.9
Uniprot ID	E6LHV2
Expiration Date	6 months from date of receipt.