

## Product Datasheet

### Cas12 (LbCpf1), Active (orb1742856)

<b>Description</b>	Cas12 (LbCpf1), Active
<b>Preservatives</b>	Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 1mM DTT, 10% glycerol.. Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 1mM DTT, 10% glycerol.
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	0.1 µg/µl
<b>Storage</b>	Store at -70°C.
<b>Tag</b>	GST
<b>Note</b>	For research use only
<b>Source</b>	Lachnospiraceae
<b>MW</b>	144.8
<b>Uniprot ID</b>	<a href="#">A0A5S8WF58</a>
<b>Expiration Date</b>	6 months from date of receipt.