

Product Datasheet

AsCas12a, Active (orb1742857)

Description	AsCas12a, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT and 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT and 10% glycerol.
Form/Appearance	liquid
Concentration	0.1 µg/µl
Storage	Store at -70°C.
Tag	His
Note	For research use only
Source	Acidaminococcus sp
MW	152.5
Uniprot ID	U2UMQ6
Expiration Date	6 months from date of receipt.