

Product Datasheet

Anti 2019-nCoV N-Protein hIgG Antibody, HRP conjugated (orb1742891)

Description	Anti 2019-nCoV N-Protein hIgG Antibody, HRP conjugated
Conjugation	Unconjugated
Preservatives	1X PBS + 0.02% sodium azide. 1X PBS + 50% glycerol + HRP stabilizer
Form/Appearance	liquid
Concentration	0.2 µg/µ
Storage	Store at -20°C. Stable for 1 year at -20°C
Note	For research use only
Source	Recombinant Human Monoclonal IgG Antibody
Expiration Date	12 months from date of receipt.