

Product Datasheet

MERS Spike Protein RBD (orb1742911)

Description	MERS Spike Protein RBD
Preservatives	Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 90mM glycine.. Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 90mM glycine.
Form/Appearance	liquid
Storage	Stored at -70°C. Stable for 1 year at -70°C.
Tag	Fc
Note	For research use only
Source	MERS-CoV (virus)
MW	50.1
Expiration Date	6 months from date of receipt.