

## Product Datasheet

# 2019-nCoV NSP14 Methyltransferase, Active (orb1742974)

<b>Description</b>	2019-nCoV NSP14 Methyltransferase, Active
<b>Preservatives</b>	Recombinant protein stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.25mM DTT, 25% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.0, 300mM NaCl, 150mM imidazole, 0.25mM DTT, 25% glycerol.
<b>Form/Appearance</b>	liquid
<b>Storage</b>	Stored at -70°C. Stable for 1 year at -70°C.
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Source</b>	2019-nCoV (virus)
<b>MW</b>	61.2
<b>Uniprot ID</b>	<a href="#">P0DTD1</a>
<b>Expiration Date</b>	6 months from date of receipt.