

Product Datasheet

**Mouse TNFR2 / CD120b / TNFRSF1B Protein, Fc Tag
(MALS verified) (orb1742992)**

Description

TNF RI is also known as the p60

p75 or p80 TNFR) are two distinct type...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TNFR2

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.8 kDa. The protein migrates as 65-70 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse TNFR2 Protein, Fc Tag (orb1742992) is expressed from human 293 cells (HEK293). It contains AA Val 23 - Gly 258 (Accession # P25119-1).

MW

51.8 kDa

Uniprot ID

P25119

Expiration Date

6 months from date of receipt.