
Product Datasheet

**Cynomolgus Apolipoprotein E Protein, His Tag
(MALS & SPR verified) (orb1742998)**

Description

Cynomolgus Apolipoprotein E

from human 293 cells (HEK293).
It contains...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Apolipoprotein E

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 35.9 kDa. The protein migrates as 35-40 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus Apolipoprotein E Protein, His Tag (orb1742998) is expressed from human 293 cells (HEK293). It contains AA Lys 19 - His 317 (Accession # P10517).

MW

35.9 kDa

Uniprot ID

P10517

Expiration Date

6 months from date of receipt.