

---

## Product Datasheet

### Cynomolgus ASGR1 Protein, His Tag (orb1743023)

**Description**

Cynomolgus ASGR1 Protein, His Tag is

(HEK293). It contains AA Ser 61 ...

<b>Reactivity</b>	Monkey
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	ASGR1
<b>Preservatives</b>	Powder
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
<b>Tag</b>	N-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 28.3 kDa. The protein migrates as 37-42 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Purity</b>	90%
<b>Source</b>	Cynomolgus ASGR1 Protein, His Tag (orb1743023) is expressed from human 293 cells (HEK293). It contains AA Ser 61 - Leu 291 (Accession # XP_005582755.1).
<b>MW</b>	28.3 kDa
<b>Uniprot ID</b>	XP_005582755.1
<b>Expiration Date</b>	6 months from date of receipt.