
Product Datasheet

Cynomolgus Flt-3 / Flk-2 Protein, His Tag (MALS verified) (orb1743025)

Description

Cynomolgus Flt-3 Protein, His Tag

(HEK293). It contains AA Gln 28 ...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Flt-3

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 59.9 kDa. The protein migrates as > 85 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus Flt-3 Protein, His Tag (orb1743025) is expressed from human 293 cells (HEK293). It contains AA Gln 28 - Ser 541 (Accession # A0A2K5VQ48).

MW

59.9 kDa

Uniprot ID[A0A2K5VQ48](https://www.uniprot.org/uniprot/A0A2K5VQ48)**Expiration Date**

6 months from date of receipt.