

---

## Product Datasheet

**Cynomolgus KLRG1 Protein, Fc Tag (SPR verified)  
(orb1743035)**

**Description**

Cynomolgus KLRG1 Protein, Fc Tag

(HEK293). It contains AA Leu 60 -...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

KLRG1

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

N-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 41.3 kDa. The protein migrates as 43-53 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Cynomolgus KLRG1 Protein, Fc Tag (orb1743035) is expressed from human 293 cells (HEK293). It contains AA Leu 60 - Pro 189 (Accession # A0A2K5WAP9).

**MW**

41.3 kDa

**Uniprot ID**[A0A2K5WAP9](#)**Expiration Date**

6 months from date of receipt.