
Product Datasheet

**Cynomolgus KLRG1 Protein, His Tag (SPR verified)
(orb1743036)**

Description

Cynomolgus KLRG1 Protein, His Tag

(HEK293). It contains AA Leu 60 ...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

KLRG1

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 16.8 kDa. The protein migrates as 18-21 kDa and 22-30 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus KLRG1 Protein, His Tag (orb1743036) is expressed from human 293 cells (HEK293). It contains AA Leu 60 - Pro 189 (Accession # A0A2K5WAP9).

MW

16.8 kDa

Uniprot ID

[A0A2K5WAP9](#)

Expiration Date

6 months from date of receipt.