

---

## Product Datasheet

### Human CCL11 Protein, His Tag (orb1743039)

**Description**

Human CCL11 Protein, His Tag is

AA Gly 24 - Pro 97 (Access...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CCL11

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag**

N-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 10.1 kDa. The protein migrates as 12 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

**Purity**

90%

**Source**

Human CCL11 Protein, His Tag (orb1743039) is expressed from E. coli cells. It contains AA Gly 24 - Pro 97 (Accession # P51671).

**MW**

10.1 kDa

**Uniprot ID**
**P51671**
**Expiration Date**

6 months from date of receipt.