

---

## Product Datasheet

**Human Beta-NGF Protein, premium grade  
(orb1743045)**

**Description**

Human Beta-NGF Protein, premium

cells (HEK293). It contains AA Ser...

**Reactivity**

Human

**Endotoxins**

0.01 EU per µg

**Conjugation**

Unconjugated

**Target**

Beta-NGF

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

Native

**Note**

For research use only

**Application notes**

This protein carries no "tag". The protein has a calculated MW of 13.5 kDa. The protein migrates as 14-15 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Human Beta-NGF Protein, premium grade (orb1743045) is expressed from human 293 cells (HEK293). It contains AA Ser 122 - Ala 241 (Accession # P01138).

**MW**

13.5 kDa

**Expiration Date**

6 months from date of receipt.