
Product Datasheet

**Human KRAS Protein, His Tag (active enzyme,
HPLC verified) (orb1743046)**

Description

Human KRAS Protein, His Tag is

 contains AA Thr 2 - Cys 185
(Accessi...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

KRAS

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 23.0 kDa. The protein migrates as 25-27 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

95%

Source

Human KRAS Protein, His Tag (orb1743046) is expressed from E. coli cells. It contains AA Thr 2 - Cys 185 (Accession # P01116-2).

MW

23.0 kDa

Uniprot ID
P01116
Expiration Date

6 months from date of receipt.