

Product Datasheet

Human CTLA-4 / CD152 Protein, Twin-Strep Tag (orb1743052)

Description

Human CTLA-4 Protein, Twin-Strep

cells (HEK293). It contains AA Ala ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CTLA-4

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Twin-Strep

Note

For research use only

Application notes

This protein carries a twin strep tag at the C-terminus. The protein has a calculated MW of 16.5 kDa. The protein migrates as 20-30 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_005205.2

Purity

85%

Source

Human CTLA-4, Strep Tag (orb1743052) is expressed from human 293 cells (HEK293). It contains AA Ala 37 - Phe 162 (Accession # P16410-1).

MW

16.5 kDa

Uniprot ID
P16410
NCBI
NP_005205.2
Expiration Date

6 months from date of receipt.