

## Product Datasheet

# Human TLR2 full length protein-synthetic nanodisc (orb1743131)

<b>Description</b>	Human TLR2 full length protein-synthetic nanodisc
<b>Reactivity</b>	Human
<b>Target</b>	TLR2
<b>Form/Appearance</b>	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Note</b>	For research use only
<b>Source</b>	Mammalian
<b>MW</b>	The human full length TLR2 protein has a MW of 89.8 kDa
<b>Uniprot ID</b>	<b>O60603</b>
<b>Expiration Date</b>	6 months from date of receipt.