

## Product Datasheet

### Human FSTL1 Protein, His Tag (orb1743266)

<b>Description</b>	Recombinant Human FSTL1 Protein with C-terminal 6×His tag
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	FSTL1
<b>Form/Appearance</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Tag</b>	C-6×His Tag
<b>Note</b>	For research use only
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>Source</b>	Mammalian
<b>MW</b>	The protein has a predicted molecular mass of 33.5 kDa after removal of the signal peptide. The apparent molecular mass of FSTL1-His is approximately 35-55 kDa due to glycosylation.
<b>Uniprot ID</b>	<a href="#">Q12841</a>
<b>Expiration Date</b>	6 months from date of receipt.