

Product Datasheet

Human IL20RA Protein, hFc Tag (orb1743311)

Description	Recombinant Human IL20RA Protein with C-terminal human Fc tag
Reactivity	Human
Conjugation	Unconjugated
Target	IL20RA
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM NaAc-HAc, pH 4.5.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-Human Fc Tag
Note	For research use only
Purity	Greater than 98% as determined by reducing SDS-PAGE.
Source	E.coli
MW	12.42 KDa
Uniprot ID	P01127
Expiration Date	6 months from date of receipt.