

## Product Datasheet

### ITGAV Antibody (orb1743393)

<b>Description</b>	Anti-ITGAV(abituzumab biosimilar) mAb
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Target</b>	ITGAV
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Note</b>	For research use only
<b>Application notes</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Isotype</b>	IgG2
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P06756</b>
<b>Dilution Range</b>	Flow Cyt 1:100
<b>Expiration Date</b>	12 months from date of receipt.