

## Product Datasheet

### SOS1 Antibody (orb1743691)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | SOS1 Antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA, FC   |
| <b>Immunogen</b>           | A synthesized peptide derived from Human SOS1   |
| <b>Target</b>              | SOS1  |
| <b>Preservatives</b>       | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                     |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | Rabbit IgG  |
| <b>Clonality</b>           | Monoclonal  |
| <b>Uniprot ID</b>          | <a href="#">Q07889</a>  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |