

Product Datasheet

CDH17 Antibody (orb1743741)

Description	CDH17 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Human CDH17 protein
Target	CDH17
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	hIgG1
Clonality	Monoclonal
Uniprot ID	Q12864
Expiration Date	12 months from date of receipt.