



Product Datasheet

Anti-NUDT6 Antibody (orb1743769)


Description	Anti-NUDT6 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E.coli-derived human NUDT6 recombinant protein (Position: D45-K303).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Mouse Direct ELISA, 0.1-0.5 µg/ml, Human. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	P53370
Expiration Date	12 months from date of receipt.



WB analysis using anti-NUDT6 antibody.1:...



IHC analysis of NUDT6 using anti-NUDT6 a...



IHC analysis of NUDT6 using anti-NUDT6 a...