

Product Datasheet

Anti-Fruit fly Rim Antibody (orb1743860)

Description	Biorbyt Anti-Rim Antibody. This antibody reacts with Fruit fly.
Species/Host	Rabbit
Reactivity	Drosophila
Conjugation	Unconjugated
Form/Appearance	Liquid
Concentration	500 µg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml
Isotype	Rabbit IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.