

Product Datasheet

Genome polyprotein Antibody (orb1743998)

Description	Genome polyprotein Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic Peptide (979-992aa)
Target	Genome polyprotein
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9WJ28
Expiration Date	12 months from date of receipt.