

Product Datasheet

Early E3 18.5 kDa glycoprotein Antibody (orb1744002)

Description	Early E3 18.5 kDa glycoprotein Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human adenovirus C serotype 5 Early E3 18.5 kDa glycoprotein protein (18-124aa)
Target	Early E3 18.5 kDa glycoprotein
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04494
Expiration Date	12 months from date of receipt.