

Product Datasheet

Tipranavir (orb1745235)

Description	Tipranavir
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{31}H_{33}F_3N_2O_5S$
Purity	99.66%
MW	602.66
CAS Number	174484-41-4
SMILES	<chem>C(CC1=CC=CC=C1)[C@]2(CCC)CC(O)=C([C@H](CC)C3=CC(NS(=O)(=O)C4=CC=C(C(F)(F)F)C=N4)=CC=C3)C(=O)O2</chem>
Expiration Date	12 months from date of receipt.