

Product Datasheet

Endomorphin 1 and 2 Antibody (orb176241)

Description	Endomorphin 1 and 2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Endomorphin-1 conjugated to BSA.
Form/Appearance	lyophilized
Storage	The product can be stored as supplied for up to 12 months at 2°C-4°C. After reconstitution, aliquot and store at -20°C for higher stability or at 4°C with an appropriate antibacterial agent. Avoid freeze thaw-cycles.
Note	For research use only
Isotype	Rabbit IgG
Dilution Range	Immunohistochemistry: 5-10 ug/ml Immunofluorescence: 5-15 ug/ml
Expiration Date	12 months from date of receipt.