

## Product Datasheet

### USP20 Immunizing peptide (orb17768)

<b>Description</b>	USP20 Immunizing peptide
<b>Preservatives</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Storage</b>	Shipped at ambient temperature, store at - 20°C
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb17768 (USP20 antibody)
<b>Protein Sequence</b>	HGEQKIEAETRAV
<b>Expiration Date</b>	12 months from date of receipt.