



Product Datasheet

Fibrinogen beta antibody (orb178377)


Description	Goat polyclonal antibody to Fibrinogen beta
Species/Host	Goat
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	fibrinogen beta chain (aa177-190)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in Human Liver lysates (calculated MW of 55.9kDa according to NP_005132.2). Recommended concentration: 0.01-0.03µg/ml. Primary incubation 1 hour at room temperature. Preliminary testing was unsuccessful on Pig Liver for this par. IHC: Paraffin embedded Human Placenta and Kidney. Recommended concentration: 5µg/ml.
Protein Sequence	EKHQLYIDETVNSN
Clonality	Polyclonal
MW	55.9; 50.0
Entrez	2244
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Pl...



Immunohistochemical staining of Human Ki...