



## Product Datasheet

### Taq Polymerase Antibody (orb1784561)

<b>Description</b>	Taq Polymerase Antibody
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HotStart PCR
<b>Immunogen</b>	Taq-Polymerase
<b>Target</b>	Taq Polymerase
<b>Preservatives</b>	Phosphate Buffered Saline (PBS) with 0.09% Sodium Azide, rAlbumin-Free
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1000 µg/ml
<b>Storage</b>	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.



Anti-Taq [TP7] decreases primer dimer fo...



Anti-Taq [TP7] inhibits Taq polymerizati...