



Product Datasheet

Taq Polymerase Antibody (orb1784562)

Description	Taq Polymerase Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	HotStart PCR
Immunogen	Taq-Polymerase
Target	Taq Polymerase
Preservatives	Phosphate Buffered Saline (PBS) with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Anti-Taq [TP7] decreases primer dimer fo...



Anti-Taq [TP7] inhibits Taq polymerizati...