



Product Datasheet

PVRL1/NECTIN1 Antibody (orb1784566)


| | |
|----------------------------|--|
| Description | PVRL1/NECTIN1 Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | FC, IP, WB |
| Immunogen | Residues 31-331 (ECD) |
| Target | PVRL1/NECTIN1 |
| Preservatives | Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free |
| Form/Appearance | Whole IgG |
| Concentration | 1000 µg/ml |
| Storage | Shelf life: 1 year from date of receipt. Storage: 2 - 8°C |
| Note | For research use only |
| Application notes | Format: Whole IgG |
| Clonality | Recombinant |
| Uniprot ID | Q15223 |
| Dilution Range | WB: 1:1000, IP: 6 µl/mg lysate |
| Expiration Date | 12 months from date of receipt. |



Detection of human PVRL1/NECTIN1 by west...



Detection of human PVRL1/NECTIN1 (shaded...



Detection of human PVRL1/NECTIN1 by west...