



## Product Datasheet

### PVRL1/NECTIN1 Antibody (orb1784567)


<b>Description</b>	PVRL1/NECTIN1 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IP, WB
<b>Immunogen</b>	Residues 31-331 (ECD)
<b>Target</b>	PVRL1/NECTIN1
<b>Preservatives</b>	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
<b>Form/Appearance</b>	Whole IgG
<b>Concentration</b>	1000 µg/ml
<b>Storage</b>	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Format: Whole IgG
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q15223</a>
<b>Dilution Range</b>	WB: 1:1000, IP: 6 µl/mg lysate
<b>Expiration Date</b>	12 months from date of receipt.



Detection of human PVRL1/NECTIN1 by west...



Detection of human PVRL1/NECTIN1 (shaded...



Detection of human PVRL1/NECTIN1 by west...