


## Product Datasheet

# Porcine CTLA4(Cytotoxic T-lymphocyte protein 4) ELISA Kit (orb1784628)

<b>Description</b>	Porcine CTLA4(Cytotoxic T-lymphocyte protein 4) ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	CTLA4
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<a href="#">Q9MYX7</a>
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Porcine  
CTLA4(Cytotoxic  
T-lymphocyte  
pro...