


## Product Datasheet

# Human anti-IL18BP(Interleukin 18 Binding Protein)antibody Elisa kit (orb1784684)

<b>Description</b>	Human anti-IL18BP(Interleukin 18 Binding Protein)antibody Elisa kit
<b>Reactivity</b>	Human
<b>Range</b>	3.125-200ng/ml
<b>Target</b>	anti-IL18BP
<b>Storage</b>	2 - 8 °C for 12 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum,Plasma,Cell Culture Supernatant,Cell or tissue lysate,Other liquid samples
<b>Uniprot ID</b>	<b>O95998</b>
<b>Sensitivity</b>	1.875ng/ml
<b>Expiration Date</b>	Please enquire.



Human anti-IL18BP(Interleukin 18 Binding Protein)antibody Elisa kit

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)