


Product Datasheet

HUMAN TRADD(TNFRSF1A-associated via death domain) ELISA Kit (orb1784689)

Description	HUMAN TRADD(TNFRSF1A-associated via death domain) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	TRADD
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15628
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



HUMAN
TRADD(TNFRSF1A-
associated via
deat...