


## Product Datasheet

# Canine IL-1Ra(Interleukin 1 Receptor Antagonist) ELISA Kit (orb1784719)

<b>Description</b>	Canine IL-1Ra(Interleukin 1 Receptor Antagonist) ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	IL-1Ra
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9BEH0</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Canine IL-1Ra(Interleukin 1 Receptor Ant...